



Ground Coffee Excel Spreadsheet Rotap Methodology

TEST #	1		2		3	
DATE:	10-Feb-03		10-Feb-03		10-Feb-03	
N O T E S	Coarse Drip Grind Sample		Fine Drip Grind Sample		Filter Fine Grind Sample (Also illustrated by Test #4)	
MESH (Tyler)	Weight (grams)	% Retained	Weight (grams)	% Retained	Weight (grams)	% Retained
10	0.6	0.6%	0.0	0.0%	0.0	0.0%
14	5.1	4.4%	2.3	2.2%	0.0	0.0%
20	34.6	28.9%	28.4	24.9%	15.2	15.0%
28	77.6	42.2%	67.0	36.8%	39.2	23.8%
PAN	102.0	23.9%	105.0	36.2%	101.2	61.3%

TEST #	4		5		6	
DATE:	10-Feb-03		10-Feb-03		10-Feb-03	
N O T E S	Filter Fine Grind Sample (Also illustrated by Test #3)		Espresso Grind Sample			
MESH (Tyler)	Weight (grams)	% Retained	Weight (grams)	% Retained	Weight (grams)	% Retained
20	15.2	15.0%	0.0	0.0%		
28	39.2	23.8%	3.2	2.9%		
35	66.4	26.9%	33.8	27.8%		
48	87.9	21.3%	75.1	37.6%		
PAN	101.2	13.1%	110.0	31.7%		